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Minsalud

INSTITUTO
NACIONAL DE
SALUD

INSTITUTO
NACIONAL DE
SALUD

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S. pneumoniae ATCC 49619


Table 2G. Zone Diameter and Minimal Inhibitory Concentration Interpretive Standards for *Streptococcus pneumoniae***Testing Conditions**

Medium: Disk diffusion: Mueller-Hinton agar (MHA) with 5% sheep blood
 Broth dilution: cation-adjusted Mueller-Hinton broth with lysed horse blood (LHB) (2.5% to 5% v/v)
 (see M07-A10 for instructions for preparation of LHB)
 Agar dilution: MHA with sheep blood (5% v/v); recent studies using the agar dilution method have not been performed and reviewed by the subcommittee.

Inoculum: Direct colony suspension, equivalent to a 0.5 McFarland standard, prepared using colonies from an overnight (18- to 20-hour) sheep blood agar plate

Incubation: 35°C ± 2°C
 Disk diffusion: 5% CO₂; 20 to 24 hours
 Dilution methods: ambient air; 20 to 24 hours (CO₂ if necessary for growth with agar dilution)

Routine QC Recommendations (See Tables 4B and 5B for acceptable QC ranges.)


Streptococcus pneumoniae ATCC® 49619

Disk diffusion: deterioration of oxacillin disk content is best assessed with *Staphylococcus aureus* ATCC® 25923, with an acceptable range of 18–24 mm on unsupplemented Mueller-Hinton agar (MHA).

When a commercial test system is used for susceptibility testing, refer to the manufacturer's instructions for QC test recommendations and QC ranges.

Antibiótico	Interpretación	Diagnóstico clínico	
		Meningitis	No meningitis
Penicilina	Sensible	$\leq 0,06 \mu\text{g/ml}$	$\leq 2,0 \mu\text{g/ml}$
	Intermedia	-	$= 4,0 \mu\text{g/ml}$
	Resistente	$\geq 0,12 \mu\text{g/ml}$	$\geq 8 \mu\text{g/ml}$
Ceftriaxona	Sensible	$\leq 0,5 \mu\text{g/ml}$	$\leq 1,0 \mu\text{g/ml}$
	Intermedia	$= 1,0 \mu\text{g/ml}$	$= 2,0 \mu\text{g/ml}$
	Resistente	$\geq 2 \mu\text{g/ml}$	$\geq 4 \mu\text{g/ml}$

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Testing Conditions	
Medium:	Disk diffusion: <i>Haemophilus</i> Test Medium (HTM) Broth dilution: HTM broth
Inoculum:	Direct colony suspension, equivalent to a 0.5 McFarland standard prepared using colonies from an overnight (preferably 20- to 24-hour) chocolate agar plate [see comment (2)]
Incubation:	35°C ± 2°C; Disk diffusion: 5% CO ₂ ; 16 to 18 hours Broth dilution: ambient air; 20 to 24 hours

Routine QC Recommendations (See Tables 4A, 4B, 5A, and 5B for acceptable QC ranges.)

Haemophilus influenzae ATCC® 49247
Haemophilus influenzae ATCC® 49766

Use either *Haemophilus influenzae* ATCC® 49247 or *Haemophilus influenzae* ATCC® 49766 or both of these strains, based on the antimicrobial agents to be tested. Neither strain has QC ranges for all agents that might be tested against *H. influenzae* or *H. parainfluenzae*.

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96

©Clinical and Laboratory Standards Institute. All rights reserved.

Table 21
Neisseria meningitidis
M02 and M07

Table 21. Zone Diameter and MIC Breakpoints for *Neisseria meningitidis*

Testing Conditions	
Medium:	Disk diffusion: MHA with 5% sheep blood Broth microdilution: CAMHB supplemented with LHB (2.5% to 5% v/v) (see M07 ¹ for preparation of LHB) Agar dilution: MHA supplemented with sheep blood (5% v/v)
Inoculum:	Colony suspension from 20–24 hours growth from chocolate agar incubated at 35°C; 5% CO ₂ ; equivalent to a 0.5 McFarland standard. Colonies grown on sheep blood agar may be used for inoculum preparation. However, the 0.5 McFarland suspension obtained from sheep blood agar will contain approximately 50% fewer CFU/mL. This must be considered when preparing the final dilution before panel inoculation, as guided by colony counts.
Incubation:	35°C ± 2°C; 5% CO ₂ ; 20–24 hours
Routine QC Recommendations (See Tables 4A-1, 4B, 5A-1, and 5B for acceptable QC ranges.)	
<i>Streptococcus pneumoniae</i> ATCC [®] 49619:	
Disk diffusion: incubate in 5% CO ₂ .	
Broth microdilution: incubate in ambient air or CO ₂ (except azithromycin QC tests that must be incubated in ambient air).	
<i>E. coli</i> ATCC [®] 25922	
Disk diffusion, broth microdilution or agar dilution for ciprofloxacin, nalidixic acid, minocycline, and sulfisoxazole: incubate in ambient air or CO ₂ .	
When a commercial test system is used for susceptibility testing, refer to the manufacturer's instructions for QC test recommendations and QC ranges.	

General Comments


Important: For complete information on safety precautions, see *Biosafety in Microbiological and Biomedical Laboratories*, 5th ed. Washington, DC: US Department of Health and Human Services; 2009. <http://www.cdc.gov/biosafety/publications/bmbl5/>. Accessed 10 December 2019.

- Recommended precautions:** Perform all AST of *N. meningitidis* in a BSC. Manipulating *N. meningitidis* outside a BSC is associated with increased risk for contracting meningococcal disease. Laboratory-acquired meningococcal disease is associated with a case fatality rate of 50%. Exposure to droplets or aerosols of *N. meningitidis* is the most likely risk for laboratory-acquired infection. Rigorous protection from droplets or aerosols is mandated when microbiological procedures (including AST) are performed on all *N. meningitidis* isolates.
- If a BSC is unavailable, manipulation of these isolates should be minimized, limited to Gram staining or serogroup identification using phenolized saline solution, while wearing a laboratory coat and gloves and working behind a full face splash shield. Use BSL-3 practices, procedures, and containment equipment for activities with a high potential for droplet or aerosol production and for activities involving production quantities or high concentrations of infectious materials. If BSL-2 or BSL-3 facilities are not available, forward isolates to a referral or public health laboratory with a minimum of BSL-2 facilities.
- Laboratorians who are exposed routinely to potential aerosols of *N. meningitidis* should consider vaccination according to the current recommendations of the Centers for Disease Control and Prevention Advisory Committee on Immunization Practices, available at <http://www.cdc.gov/vaccines/acip/index.html>.

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